

# Applications of High-Resolution MFM System with Low Moment Probe and Q-control in a Vacuum



T. Yamaoka<sup>1</sup>, K. Watanabe<sup>1</sup>, Y. Shirakawabe<sup>1</sup>, K. Chinone<sup>1</sup>, E. Saito<sup>2</sup>, M. Tanaka<sup>2</sup>, H. Miyajima<sup>2</sup>  
 (1) SII NanoTechnology Inc. 261-8507 Chiba, Japan (2) Dept. of Phys., Keio Univ. 270-2222 Yokohama, Japan

## Abstract

Magnetic force microscopy (MFM) is very useful for observing magnetic domain structures. However, due to stray fields from a MFM probe, observations of small magnetic domain structures are limited. We have developed a high-resolution MFM system that utilizes a low moment probe and a Q-controlled probe-driver, which allows the sensitive measurement in a vacuum without disturbing domain structures. Using this system, resolution finer than 20nm was achieved. Here, advantages of this MFM are demonstrated using a honeycomb nano-network, a semicircular loop and a cross-shaped pattern.

## 1. What is MFM?

**MFM image**  
 Attractive (Phase shift: Large)  
 Repulsive (Phase shift: Small)

**Multi domain** **Vortex** **Single domain**

Sample: Dr. Nakatani at NIMS

## 2. Unsolved Problems and Solution

**Conventional Probe (High Moment Probe)** vs **Low Moment Probe**

Stray fields from conventional MFM probes often disturb domain structures of a sample.

Thickness of magnetic coat: **Thinner**

**Advantages**  
 Non-disturb  
 High-Resolution  
 Easy manufacturing

## 3. Low Moment Probe<sup>[1]</sup>

CoPtCr sputtering  
 Tip(Si)  
 diameter  
 thickness  
 Cantilever  
 SEM observation

M-H curve measurement: Ishio Lab., Akita University  
 $H_c = 600 \text{ [Oe]}$


[A] t = 72(nm) d = 70nm  
 [B] t = 48(nm) d = 50nm  
 [C] t = 24(nm) d = 30nm  
 [D] t = 12(nm) d = 20nm  
 [E] Non-coat d = 10nm

Conventional Probe vs Low Moment Probe


Non-disturb  
 Small MFM sensitivity

## 4. MFM sensitivity

What is Q-factor ?



MFM sensitivity



Q-factor


$$\Delta\phi = \frac{Q}{K} \cdot \frac{\partial F}{\partial z}$$

: MFM signal (deg)  
 Q: Q-factor  
 K: Spring constant of cantilever (N/m)  
 F/z: Force gradient (N/m)

There is a range of Q-factor with the stable and sensitive.

## 5. What is Q-control?

Principle of Q-control



B. Anczykowski et al, Appl. Phys. A66, pp.885-889, 1998.  
 A.D.L. Humphris et al, Surface Science, Vol.491, No.3, pp.468-472, 2001.

Equation of motion

$$m\ddot{Z}(t) + \eta\dot{Z}(t) + kZ(t) = F_0 e^{i\omega t} + Ge^{i\frac{\omega}{2}Z(t)}$$


$$m\ddot{Z}(t) + \eta_{\text{eff}}\dot{Z}(t) + kZ(t) = F_0 e^{i\omega t}$$

Effective viscous coefficient  
 $\eta_{\text{eff}} = \eta - \frac{G}{\omega}$


Effective Q-factor  
 $Q_{\text{eff}} = \frac{m\omega}{\eta_{\text{eff}}} = \frac{m\omega}{\eta - \frac{G}{\omega}}$

## 6. Comparison of Q-control in Air / in a Vacuum

Changes of Q-curve when approaching the probe close to the sample




Distance dependence with Q-factor



Q-curve and Q-factor in air is influenced by the change of the viscous drag.

This phenomenon isn't desirable in MFM measurement.

Relation between Q-factor and MFM sensitivity



Cantilever  
 Resonance frequency 250 ~ 300kHz  
 Spring constant 40N/m  
 Low moment probe (t=24nm)

Sample  
 Perpendicular recording medium  
 200kFCI (Bit length: 127nm)

Range of stable Q  
 in air: 200 Q 2000  
 in a vacuum: 2000 Q 8000

Q-control in a vacuum is higher sensitivity and wider range of stable Q-factor<sup>[1]</sup>.

# Applications of High-Resolution MFM System with Low Moment Probe and Q-control in a Vacuum




T. Yamaoka<sup>1</sup>, K. Watanabe<sup>1</sup>, Y. Shirakawabe<sup>1</sup>, K. Chinone<sup>1</sup>, E. Saito<sup>2</sup>, M. Tanaka<sup>2</sup>, H. Miyajima<sup>2</sup>  
 (1) SII NanoTechnology Inc. 261-8507 Chiba, Japan (2) Dept. of Phys., Keio Univ. 270-2222 Yokohama, Japan

## 7. High-Resolution MFM System




Environment Controllable SPM  
 ( In air, in a vacuum and in liquid )  
 P :  $10^{-5}$  (Pa) by a turbo molecule pump  
 T : -140 ~ +800 ( ) by LN<sub>2</sub> and a heater  
 SPI4000/ E-sweep and SPA300HV,  
 SII NanoTechnology Inc.


## 8. Perpendicular Recording Media<sup>[1]</sup>




## 10. NiFe Semicircular Wire Loop<sup>[3]</sup>




### High-Resolution MFM observations



## 9. NiFe Honeycomb Nano-Network<sup>[2]</sup>



## 11. NiFe Cross-Shaped Pattern



### Summary

We have demonstrated the effectiveness of the high-resolution MFM that involves the use of a low moment probe and Q-control system in a vacuum. Standard high moment probe disturbs the domain structure of the sample. In contrast, the low moment probe allows the stable and clear observation of magnetic domain. The high-resolution MFM with low-moment probe and Q-control in a vacuum should be powerful for exploring nano-scale magnetism.

### Acknowledgment

The authors wish to thank Prof. Ishio at Akita University for useful discussion about the MFM technology and for M-H loop measurements

### Reference

- [1] T. Yamaoka, K. Watanabe, Y. Shirakawabe, and K. Chinone, "High-Sensitivity, High-Resolution Magnetic Force Microscopy System," J. Magn. Soc. Jpn., vol.27, pp.429-433, 2003.
- [2] E. Saitoh, M. Tanaka, H. Miyajima, and T. Yamaoka, "Domain-wall trapping in a ferromagnetic nano-wire network," J. Appl. Phys., vol.93, no.10, pp.7444-7446, 2003.
- [3] E. Saitoh, H. Miyajima, T. Yamaoka, and G. Tatara, "Current-induced resonance and mass determination of a single magnetic domain wall," Nature, vol.432, pp.203-206, 2004.

### High-Resolution MFM observations

